

Product datasheet for **MR204174A1V**

Mouse Zdhhc3 (NM_026917) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Zdhhc3 (NM_026917) AAV Particle
Tag:	Myc-DDK
Symbol:	Zdhhc3
Synonyms:	1110020O22Rik; 1810006O10Rik; 2210017C02Rik; DHHC-3; GODZ; Gramp1; Zfp373
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR204174 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

```

ATGATGCTTATCCCCACCCATCACTTCCGAGACATTGAGCGGAAACCAGAGTACCTCCAGCCAGAGAAGT
GTGCCCCACCTCCCTTCCCTGGTCTGCGGGAGCCATGTGGTTTATCCGAGATGGCTGTGGCATTGCTTG
TGCCATTGTACCTGGTTTCTGGTCTCTATGCGGAGTTTGTAGTCCTCTTTGTCATGCTGGTTCATCC
CGAGACTACGCGTACAGCATCATCAACGGAATTGTGTTCAACCTGCTGGCCTTCTTGGCCTGGCCTCCC
ACTGCCGGGCCATGCTGACGGACCCCGGGCAGTGCCCAAAGGAAATGCCACTAAAGAGTTCATCGAGAG
CCTTCAGCTGAAGCCTGGGCAGGTGGTGTACAAGTGTCCTCAAGTGTGCAGCATCAAGCCCGACCGGGCA
CACCAGTGCAGTGTTTGTAAAGCGGTGCATTGCAAGATGGATCACCAGTGTCTTGGGTCAACAAGTGTG
TCGGCGAGAACAACCAGAAGTACTTTGTCTATTACAAATGTACATAGCTCTCATTTCTTGCACGCCCT
CATCATGGTGGGATTCCACTTCTGCTGCTTTGAAGAAGACTGGACAAAGTGCAGTCTCTTCTCGCCG
CCCACCACAGTCATCCTGCTCATCTGCTGTGCTTTGAGGCCCTGCTCTTCTCATTTTCACATCAGTGA
TGTTTGGGACCAAGTACACTCCATCTGCACAGATGAGACGGGCATAGAACAATTGAAAAAGGAAGAGAG
AAGATGGGCTAAAAAACAATAAGATGAACATGAAAGCCGTGTTGGCCACCCTTCTCTTAGGCTGG
GCCAGCCCTTTGCCACACCAGACCAAGGAAGGCAGACCCGTACCAGTATGTGGTC

```

Protein Sequence: >MR204174 protein sequence
 Red=Cloning site Green=Tags(s)

```

MMLIPTHHFRDIERKPEYLQPEKCAPPPFPGPAGAMWFIRDGCGIACAIVTWFLVLYAEFVVLVFMVLVPS
RDYAYSIIINGIVFNLLAFLALASHCRAMLTDPGAVPKGNATKEFIESLQLKPGQVYKCPKCCSIKPDRA
HHCSVKRCIRKMDHHCWPVNNVCVGENNQKYFVLFMYIALISLHALIMVGFHFLHCFEEDWTKCSSFSP
PTTVILLILLCFEALLFLIFTSVMFGTQVHSICTDETGIEQLKKEERRWAKKTKWMNMKAVFGHPFSLGW
ASPFATPDQ GKADPYQYVV

```


[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_026917
ORF Size:	840 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_026917.4
RefSeq Size:	6987 bp
RefSeq ORF:	900 bp
Locus ID:	69035
UniProt ID:	Q8R173
Cytogenetics:	9 F4
MW:	34 kDa