

## Product datasheet for **MR204180L3V**

### **Slc25a35 (NM\_028048) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Slc25a35 (NM_028048) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc25a35
Synonyms:	1810012H11Rik
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
ACCN:	NM_028048
ORF Size:	903 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204180).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_028048.1</a>
RefSeq Size:	1827 bp
RefSeq ORF:	903 bp
Locus ID:	71998
UniProt ID:	<a href="#">Q5SWT3</a>



[View online »](#)