

## Product datasheet for **MR204183L3V**

### Fam212b (NM\_175398) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Fam212b (NM_175398) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fam212b
Synonyms:	6530418L21Rik; inka2
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
ACCN:	NM_175398
ORF Size:	903 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204183).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_175398.3</a>
RefSeq Size:	4701 bp
RefSeq ORF:	903 bp
Locus ID:	109050
UniProt ID:	<a href="#">Q80VY2</a>
Gene Summary:	Inhibitor of the serine/threonine-protein kinase PAK4. Acts by binding PAK4 in a substrate-like manner, inhibiting the protein kinase activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)