

## Product datasheet for MR204205

### Pou6f1 (NM\_010127) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Pou6f1 (NM\_010127) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Pou6f1  
**Synonyms:** 2310038G18Rik; C79593; cns-1; Emb  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR204205 representing NM\_010127  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGCCCGGGATCAGCAGTCAGATCCTCACGAATGCTCAGGGACAGGTTATTGGAGCACTCCGTGGGTAG  
 TGAAGTCACTAGCGTGGCCACACCAGCACCGGCACAGAGCTGCAGGTCCAAGCCGTGACTCCCCAGT  
 CTTGTTGAATGCCAGGGCCAGGTGATCGCAACCCTAGCCAGCAGCCCCCTGCCTCAGCCTGTGGCTGTC  
 AGGAAGCCAAACACACCCGAGTCCCCTGCTAAGAGTGAGGTGCAGCCTATCCAGCCGACACAAGCCGTGC  
 CCCAGCCTGCAGTAATCCTCACCAGCCCAACGCCAGCGCTCAAGCCGTGACTGCAACTCCCATCCCAAT  
 CACCTGCTCAGAGACCCCAACCGTCAGTCAGTTGGTATCAAAGCCGCACACCCCAAGTCTGGATGAGGAC  
 GGGATCAACTTAGAAGAGATCCGGGAGTTTGCTAAGAATTTAAGATCCGGCGGCTCTCCCTGGGTCTGA  
 CACAGACCCAGGTGGGCCAGGCTTTGACGGCGACAGAAGGGCCGGCTACAGCCAATCAGCCATCTGCAG  
 GTTTGAGAAATTGGACATCACACCAAGAGCGCCAGAGCTGAAGCCGGTTTTGGAAAAGTGGTTGATG  
 GAGGCAGAGCTCCGCAACCAGGAAGGCCAGCAGAACCTGATGGAGTTCGTGGCGGCGAGCCCTCCAAGA  
 AACGCAAGCGGCGCACTTCCTTCACACCGCAGGCCATAGAGGCTCTCAATGCCTACTTTGAGAAAAACCC  
 CCTGCCACCGGCCAGGAGATCACGGAGATCGCCAAGGAGCTCAACTACGACCGGAGGTGGTGAAGGTC  
 TGTTCTGTAATCGACGCCAGACACTGAAGAACCAGCAAGCTGAACGTCTTTCAGATCCCG

**ACGCGT**ACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR204205 representing NM\_010127  
Red=Cloning site Green=Tags(s)


MPGISSQILTNAQGQVIGALPWVNSASVATPAPAQSLQVAVTPQLLLNAQGQVIATLASSPLPQPVAV  
 RKPNTPEPAKSEVQPIQPTQAVPQPAVILTSPTPALKPSAATPITCSETPTVSQLVSKPHTPSLDED  
 GINLEEIREFAKNFKIRRLSLGLTQTQVGQAL TATEGPAYSSQSAICRFEKLDITPKSAQKLPVLEKWL  
 MAELRNQEGQQLMEFVGGPEPSKRRKRRTSFTPQAIEALNAYFEKNPLPTGQEITEIAKELNYDREVV  
 WFCNRRQTLKNTSKLNVFQIP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9058\\_d04.zip](https://cdn.origene.com/chromatograms/mm9058_d04.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**




**ACCN:** NM\_010127

**ORF Size:** 903 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_010127.3, NP_034257.1</u>
<b>RefSeq Size:</b>	4575 bp
<b>RefSeq ORF:</b>	906 bp
<b>Locus ID:</b>	19009
<b>UniProt ID:</b>	<u>Q07916</u>
<b>Cytogenetics:</b>	15 56.35 cM
<b>MW:</b>	32.9 kDa
<b>Gene Summary:</b>	Transcription factor that binds preferentially to a variant of the octamer motif (5'-ATGATAAT-3').[UniProtKB/Swiss-Prot Function]

**Product images:**


Circular map for MR204205