

Product datasheet for **MR204214L4V**

Lenti ORF particles, MGC:38152 (GFP-tagged) - Mouse cDNA clone MGC:38152 IMAGE:5321417, 200ul, >10⁷ TU/mL (Caprin2) (BC025636) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lenti ORF particles, MGC:38152 (GFP-tagged) - Mouse cDNA clone MGC:38152 IMAGE:5321417, 200ul, >10 ⁷ TU/mL (Caprin2) (BC025636) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Caprin2
Synonyms:	C1qdc1; Eeg1; rng140
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC025636
ORF Size:	903 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204214).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC025636.1
RefSeq Size:	3667 bp
RefSeq ORF:	905 bp
Locus ID:	232560
Cytogenetics:	6 G3



[View online »](#)

Gene Summary:

This gene encodes a member of the caprin family. The encoded protein may function as an RNA-binding protein that induces the formation of RNA granules and plays an important role in brain function. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2014]