

Product datasheet for **MR204223L3V**

Dhrs3 (NM_001172424) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dhrs3 (NM_001172424) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dhrs3
Synonyms:	retSDR1; Rsdr1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001172424
ORF Size:	909 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204223).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001172424.1 , NP_001165895.1
RefSeq Size:	1367 bp
RefSeq ORF:	831 bp
Locus ID:	20148
Cytogenetics:	4 E1
Gene Summary:	Catalyzes the reduction of all-trans-retinal to all-trans-retinol in the presence of NADPH. [UniProtKB/Swiss-Prot Function]



[View online »](#)