

Product datasheet for **MR204239L4V**

Rpp40 (NM_145938) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rpp40 (NM_145938) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rpp40
Synonyms:	D8Bwg1265e; Rnasep1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145938
ORF Size:	912 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204239).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145938.2
RefSeq Size:	2736 bp
RefSeq ORF:	1092 bp
Locus ID:	208366
UniProt ID:	Q8R1F9
Cytogenetics:	13 A3.3
Gene Summary:	Component of ribonuclease P, a ribonucleoprotein complex that generates mature tRNA molecules by cleaving their 5'-ends. Also a component of the MRP ribonuclease complex, which cleaves pre-rRNA sequences.[UniProtKB/Swiss-Prot Function]



[View online »](#)