

## Product datasheet for **MR204254L4V**

### **Rdh12 (BC016204) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Rdh12 (BC016204) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rdh12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC016204
ORF Size:	912 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204254).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC016204</a> , <a href="#">AAH16204</a>
RefSeq Size:	1659 bp
RefSeq ORF:	914 bp
Locus ID:	77974
Cytogenetics:	12 C3
Gene Summary:	The protein encoded by this gene is an NADPH-dependent retinal reductase whose highest activity is toward 9-cis and all-trans-retinol. The encoded enzyme also plays a role in the metabolism of short-chain aldehydes but does not exhibit steroid dehydrogenase activity. Defects in the human gene are associated with Leber congenital amaurosis type 13, and Retinitis Pigmentosa 53. [provided by RefSeq, Sep 2015]



[View online »](#)