

Product datasheet for **MR204263L3V**

Tcrb (BC030075) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Tcrb
Synonyms:	BVI3S1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC030075
ORF Size:	912 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR204263).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	BC030075
RefSeq Size:	1125 bp
RefSeq ORF:	914 bp
Locus ID:	269846
Cytogenetics:	6 19.71 cM 6 B1

