

## Product datasheet for **MR204286L3V**

### Cnn2 (NM\_007725) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | Cnn2 (NM_007725) Mouse Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | Cnn2   |
| Synonyms:                 | AA408047; AI324678; Calpo2   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_007725  |
| ORF Size:                 | 918 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(MR204286).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_007725.2</a> , <a href="#">NP_031751.1</a>  |
| RefSeq Size:              | 1800 bp  |
| RefSeq ORF:               | 918 bp   |
| Locus ID:                 | 12798  |
| UniProt ID:               | <a href="#">Q08093</a>   |
| Cytogenetics:             | 10 C1  |



[View online »](#)

**Gene Summary:**

Thin filament-associated protein that is implicated in the regulation and modulation of smooth muscle contraction. It is capable of binding to actin, calmodulin, troponin C and tropomyosin. The interaction of calponin with actin inhibits the actomyosin Mg-ATPase activity.[UniProtKB/Swiss-Prot Function]