

Product datasheet for MR204296L3

Eif2b1 (NM_145371) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Eif2b1 (NM_145371) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Eif2b1
Synonyms:	26kD; D5Ertd406; D5Ertd406e; EIF2; EIF2BA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204296).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_145371
ORF Size:	915 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_145371.4 , NP_663346.1
RefSeq Size:	2012 bp
RefSeq ORF:	918 bp
Locus ID:	209354
UniProt ID:	Q99LC8
Cytogenetics:	5 63.67 cM
Gene Summary:	This gene encodes the alpha subunit of the eukaryotic translation initiation factor complex 2B (eIF2B). The eIF2B complex is a heterodecamer comprised of two molecules each of alpha, beta, gamma, delta and epsilon subunits. The eIF2B complex is a critical regulator of protein synthesis acting as the guanine nucleotide exchange factor for eIF2 to enable the formation of ternary complex that is required for the initiation of mRNA translation. [provided by RefSeq, Aug 2015]

Product images:



Circular map for MR204296L3