

Product datasheet for **MR204310L3V**

Armc10 (NM_026034) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Armc10 (NM_026034) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Armc10
Synonyms:	2810037C14Rik; Svh
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026034
ORF Size:	921 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204310).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026034.1
RefSeq Size:	2014 bp
RefSeq ORF:	921 bp
Locus ID:	67211
UniProt ID:	Q9D0L7
Cytogenetics:	5 A3
Gene Summary:	May play a role in cell survival and cell growth. May suppress the transcriptional activity of p53/TP53 (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)