

Product datasheet for **MR204343L3V**

Tyms (BC020139) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tyms (BC020139) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tyms
Synonyms:	TS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC020139
ORF Size:	921 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204343).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC020139 , AAH20139
RefSeq Size:	986 bp
RefSeq ORF:	923 bp
Locus ID:	22171
Cytogenetics:	5 15.81 cM



[View online »](#)

Gene Summary:

This gene encodes an enzyme that catalyzes the methylation of deoxyuridylate to deoxythymidylate using 5,10-methylenetetrahydrofolate as a cofactor. This function maintains the thymidine-5-prime monophosphate concentration critical for DNA replication and repair. The encoded enzyme is a target for cancer chemotherapeutic agents. The majority of transcripts for this gene lack a 3' UTR (PMID: 3022294, 3444407). The stop codon in these transcripts is UAA, compared to the UAG found in the genome and longer transcripts, as the polyA site is located within the stop codon (PMID: 3444407, 2157203). A related pseudogene has been identified on chromosome 10. [provided by RefSeq, Mar 2010]