

## Product datasheet for MR204355

### Ftsj1 (BC011144) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Ftsj1 (BC011144) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ftsj1
Synonyms:	Ftsjl, Sfc12, Ftsj
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR204355 representing BC011144 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGGACGGACATCTAAGGATAAACGGGATGTCTACTACCGCCTGGCCAAGGAGAATGGCTGGCGGGCCC  
GAAGCGCTTTCAAGCTCCTGCAACTGGATGAGGAATTCTCGTGCTAGAGAGTGTGAAGCGAGCGGTTGC  
AGTGTGTGCAGCCCTTTTTTTTTTTTTTTTGTGCGTCTGAAAAGCGGGGGCCACGGTTCTGGTCAG  
GTGGTGGCCGTGGACTTACAGGCTATGGCTCCACAACCAGGTGTGATACAGATACAGGGCGACATCACTC  
AGCTTTCAACTGCCAAGGAGATCATCCAGCACTTTGAAGGCTGCCCTGCTCTGTTTGTGTGTGAAGG  
GGCTCCTGATGTCACTGGTCTCCATGATGTGGATGAGTACATGCAGGCCAGCTTCTGCTAGCTGCTCTC  
AACATCGTACTCATGTCTTGAAGCTAGGGGGCTGCTTCGTGGCCAAGATATCCGAGGCCGGGATGTGA  
CCCTGCTCTACAGCCAACACTACGTATCTTCTTCCCCAGTGTGCTCTGCGCTAAGCCCAAGAGCAGCCGGAA  
CTCTAGCATTGAGGCCTTTGCTGTGTGTCAGGGTTATGACCCCTCTGAGGGCTTCATTCCAGACCTGACC  
AGACCCCTGTGAACCACTCATATGATTCAACCAGTTGGATGGTCCCACCCGTGTTATTGTGCCCTTTG  
TGGCCTGTGGGGACCTCAGCGCTACGATTCGGACCGCACTTACTCTCTGGATCCCCCATCGCACCCCC  
ATACCAAGAGGCCTGCAGGTTGAAGAAGAATGGACAGCTGGCCAAGGAACTCTCCCCAAGAATGCTCC  
ATCAACAGTGTGGACAAGTTGCCCGAGCCTCTGGCTATCCATACCCTGTAGACCCCAAGGTGGAAGACA  
ATGAAATTCATTGT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC011144.1](#)

**RefSeq Size:** 1543 bp


**RefSeq ORF:** 926 bp

**Locus ID:** 54632

**Cytogenetics:** X 3.73 cM

**MW:** 56.5 kDa

**Product images:**



Circular map for MR204355