

Product datasheet for **MR204411A1V**

Mouse Hs3st1 (NM_010474) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Hs3st1 (NM_010474) AAV Particle
Tag:	Myc-DDK
Symbol:	Hs3st1
Synonyms:	3-Ost; D5Wsu110e; Hsg3ost
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR204411 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACCTTGCTGCTCCTGGGTGCGGTGCTGCTGGTGGCCAGCCCCAGCTTGTGCATTCCCACCCGGCTG
CTCCTGGCCCGGGGCTCAAACAGCAGGAGCTTCTGAGGAAGGTGATTATTCTCCAGAGGACACCGGAGA
AGGCACAGCATCCAATGGTTCACACAGCAGCTGCCACAGACCATCATCATTGGGGTGCGAAGGTGCT
ACCCGAGCCCTGCTAGAGATGCTCAGCCTGCATCCTGATGTTGCTGCAGCTGAAAACGAGGTCCATTTCT
TTGACTGGGAGGAGCATTACAGCCAAGGCCTGGCTGGTACCTCACCCAGATGCCCTTCTCCTCCCTCA
CCAGCTCACCGTGGAGAAGACACCCGCCTATTTCACTTCGCCCAAAGTGCCTGAGAGAATCCACAGCATG
AACCCCAACCATCCGCCTGCTGCTTATCCTGAGGGACCCATCAGAGCGCGTGCTGTCCGACTACACCCAGG
TGTTGTACAACCACTTCAGAAGCACAAGCCCTATCCACCCATTGAGGACCTCCTAATGCGGGACGGTCG
GCTGAACCTGGACTACAAGGCTCTCAACCGCAGCCTGTACCATGCACACATGCTGAACCTGGCTGCGTTTT
TTCCCGTTGGGCCACATCCACATTGTGGATGGCGACCGCCTCATCAGAGACCTTTCCCTGAGATCCAGA
AGGTCGAAAGATTCTGAAGCTTTCTCCACAGATCAACGCCTCGAACTTCTACTTTAACAAAACCAAGGG
CTTCTACTGCCTGCGGGACAGTGGCAAGGACCGCTGCTTACACGAGTCCAAAGGCCGGGCGCACCCCCAG
GTGGATCCAACTACTTGATAAACTGCACGAATACTTTCATGAGCCAAATAAGAAATTTTCAAGCTCG
TGGGCAGAACATTCGACTGGCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR204411 protein sequence Red=Cloning site Green=Tags(s) MTLLLLGAVLLVAQPQLVHSHPAAPGPGGLKQQELLRKVIILPEDTGEGTASNGSTQQLPQTIIIGVRKGG TRALLEMLSLHPDVAAAENEVHFFDWEHYSQGLGWYLTQMPFSSPHQLTVEKTPAYFTSPKVPERIHS NPTIRLLLILRDPSEVLSDYTQVLYNHLQKHKPYPIEDLLMRDGRNLNDYKALNRSLYHAHMLNWLRF FPLGHIHIVDGDRLIRDPPFEIQKVERFLKLSPIQINASFYFNKTKGFYCLRDGSGKDRCLHESKGRAHPQ VDPKLLDKLHEYFHEPNKKFFKLVGRTFDWH TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_010474
ORF Size:	933 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_010474.1</u>
RefSeq Size:	1685 bp
RefSeq ORF:	936 bp
Locus ID:	15476
UniProt ID:	<u>O35310</u>
Cytogenetics:	5 21.14 cM
MW:	35.9 kDa