

## Product datasheet for **MR204411L3V**

### Hs3st1 (NM\_010474) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hs3st1 (NM_010474) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hs3st1
Synonyms:	3-Ost; D5Wsu110e; Hsg3ost
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010474
ORF Size:	936 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204411).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010474.1</a>
RefSeq Size:	1685 bp
RefSeq ORF:	936 bp
Locus ID:	15476
UniProt ID:	<a href="#">O35310</a>
Cytogenetics:	5 21.14 cM


[View online »](#)

**Gene Summary:**

Sulfotransferase that utilizes 3'-phospho-5'-adenylyl sulfate (PAPS) to catalyze the transfer of a sulfo group to position 3 of glucosamine residues in heparan. Catalyzes the rate limiting step in the biosynthesis of heparan sulfate (HSact). This modification is a crucial step in the biosynthesis of anticoagulant heparan sulfate as it completes the structure of the antithrombin pentasaccharide binding site.[UniProtKB/Swiss-Prot Function]