

## Product datasheet for MR204428L4

### Bfar (NM\_025653) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Bfar (NM_025653) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Bfar
Synonyms:	3010001A07Rik; 3110001I22Rik; AI666707; AW107665; Bar; Rnf47
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204428).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.

ACCN:	NM_025653
ORF Size:	939 bp



[View online »](#)

