

Product datasheet for **MR204431L3V**

Ctso (NM_177662) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ctso (NM_177662) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ctso
Synonyms:	A330105D01Rik; AI118514
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_177662
ORF Size:	936 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204431).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177662.2 , NP_808330.1
RefSeq Size:	3453 bp
RefSeq ORF:	939 bp
Locus ID:	229445
UniProt ID:	Q8BM88
Cytogenetics:	3 E3
Gene Summary:	This gene encodes a member of the cathepsin family of cysteine proteases that are involved in the degradation of cellular proteins. The encoded preproprotein undergoes proteolytic processing to generate a mature, functional enzyme. [provided by RefSeq, Jan 2016]



[View online »](#)