

Product datasheet for **MR204454L4V**

Qpct (BC020023) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Qpct (BC020023) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Qpct
Synonyms:	5730422A13Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC020023
ORF Size:	939 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204454).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC020023 , AAH20023
RefSeq Size:	1486 bp
RefSeq ORF:	941 bp
Locus ID:	70536
Cytogenetics:	17 E3
Gene Summary:	Responsible for the biosynthesis of pyroglutamyl peptides. Has a bias against acidic and tryptophan residues adjacent to the N-terminal glutaminyl residue and a lack of importance of chain length after the second residue (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)