

Product datasheet for **MR204480A1V**

Mouse Lgals12 (NM_019516) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Lgals12 (NM_019516) AAV Particle
Tag:	Myc-DDK
Symbol:	Lgals12
Synonyms:	AW987437; gal-12; Grip1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR204480 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCAACTGACGAACACCTGGACCCGATTCTGACAGCTTCATCCTGCAGCCGCCAGTCTTCCACCCGG
TGATTCCTTATGGCACAACAATTTTGGTGGCCTGTATGCAGGCAAGATGGTCACGCTGCAGGGTGTGGT
CCCTCTGCATGCAAGGAGGTTTCAGGTGGACTTCCAGTGTGGGTGCTGCCTGCATCCTCAGCCAGATGTT
GCCTTCCGCTTCAGCCCTCGATTCTACACTGTCAAGCCCCATGTCATCTGCAACACCCACCAAGGTGGAC
TCTGGCAAAAAGAGATCCGGTGGCCAGGGTGCCTTCAGAGAGGGGATAGCTTCCTCATTCTCTTTCT
CTTTGAGAACGAAGAGGTGAAGGTGAGTGTAAATGGCCAGCACTTTCTTCACTACCGCTACCGGCTCCCA
CTGTACGGGTAGATACCCTTGACATATCTGGTGACATCTTGGTAAAGGCTGTTGGATTCTGAACATCA
ATCCATTTGTGGAGGGTAGCAGAGAGTATCCAGTTGGATATCCCTTCTGCTGTATAGCCCCAGGCTGGA
GGTGCCCTGCTCACGTGCTCTTCTCGGGGTCTCTGGCCTGGGCAAGTCATTGTAGTTCGAGGACTGGTC
TTGAAAGAGCCGAAAGATTTTACCCTGAGCCTGAAGGATGGGACCACCCATGTTCTGTGACACTCAGGG
CTTCCTTCACAGACAGAACTGGCCTGGGTCTCCTCCTGGGACGAAAGAAGCTGATCTCGCCCCCTT
CCTCTTTCACCCCAGCGATTCTTCAGGTACTGCTTCTGTGCCAAAAGGGAGGGCTGAAGTTGGCACTC
AATGGGACAGGGCTGGGGGCCACCGCTGGACCAGAAAGCCCTGGAGCAGCTGCGGGAGCTCAGGATCA
GTGGAATGTCCACCTCTACTGTGTCCACTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR204480 protein sequence Red=Cloning site Green=Tags(s) MSTDEHLDPIPDSFILQPPVFHPVIPYGTTFGGLYAGKMVTLQGVVPLHARRFQVDFQCGCCLHPQPDV AFRFSRPFYTVKPHVICNTHQGGLWQKEIRWPGVALQRGDSFLILFLFENEVKVSVNGQHFLHYRYRLP LSRVDTLDISGDILVKAVGFLNINPFVEGSREYPVGYPFLLYSPRLEVPCSRALPRGLWPGQVIVVRGLV LKEPKDFTLSLKDGTTHVPVTLRASFTDRTLAWVSSWGRKKLISAPFLFHPQRF FEVLLLCQKGGLKLAL NGQGLGATSLDQKALEQLRELISGNVHLYCVHC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_019516
ORF Size:	942 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_019516.2, NP_062389.1</u>
RefSeq Size:	2755 bp
RefSeq ORF:	945 bp
Locus ID:	56072
UniProt ID:	<u>Q91VD1</u>
Cytogenetics:	19 A
MW:	35.5 kDa