

## Product datasheet for **MR204483L3V**

### Slc25a44 (NM\_001145876) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc25a44 (NM_001145876) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc25a44
Synonyms:	6720482A19; B430110G05Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145876
ORF Size:	942 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204483).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001145876.1</a> , <a href="#">NP_001139348.1</a>
RefSeq Size:	3546 bp
RefSeq ORF:	945 bp
Locus ID:	229517
UniProt ID:	<a href="#">Q8BGF9</a>
Cytogenetics:	3 F1



[View online »](#)

**Gene Summary:**

Mitochondrial solute transporter which transports branched-chain amino acid (BCAA; valine, leucine and isoleucine) into mitochondria in brown adipose tissue (BAT) (PubMed:31435015). BAT is involved in BCAA catabolism and actively utilizes BCAA in the mitochondria for thermogenesis (PubMed:31435015).[UniProtKB/Swiss-Prot Function]