

Product datasheet for **MR204497A1V**

Mouse Ccdc92 (NM_144819) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Ccdc92 (NM_144819) AAV Particle
Tag:	Myc-DDK
Symbol:	Ccdc92
Synonyms:	AW050178; D5Bwg0834e
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR204497 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCAGCCACAAACCTGGAGAACCAGCTGCACAGTGCACAGAAGAACCTGCTCTTCTCCAGCGGGAGC
ACGCCAGCAGCACTCAAGGGGCTGCACGCCGAGATCCGGCGGCTGCAGCAGCACTGCACAGATCTAACATA
TGAGCTGACACTCAAAAGTTTCGAACTGACAGGAGACAGCTCTCCAGAACGACGGAGCTGAAGAGGCGC
TGTGAGGAGCTGGAAGCCAGCTGAAGGCCAAAGAGGAGGAAAACCGCAAGCTGCTGCAGGAACCTGGAAC
AGAAGAACGCCGCCATCGCTGTGCTGGAGAACACAGTCCGTGAGAGAGAGAAGAAGTACCTGGAGGAGCT
GAAGGTGAAGAGCCACAACTGAGCATGCTGTGAGGTGAGTTGGAGCAACGCGCCAGCACTGTGGCCTAC
CTCACCTCGCAGCTGCATGCAGCCAAGAAGAAGTTGCTGAGCTCAAGTGGCACCTCGGATGCCAGCCCTG
CCGGCAGTCCCGCGCTGGCCAGCTACAAACCCACACCCCCAAGGACAAGCTGCCTGAGACACCCCGCCG
CCGCATGAAGAAGAGCCTGTCCGCCCACTGCACCCCGAGTTTGAGGAGGTCTACAGGTTTGAGCCGAG
AGCCGCAAACTCCTCTTCGGGAGCCTGTGGATGCCATGCCGACCCACCCATTCTGCTGGCCCGGG
AGTCTGCTGAGGTCCAGCTCAAGGAGCGGCTCTCGTCATCCCTCCCATTGCCTCAGATCGTAGCGCCAC
TGGGCAGCACAGCCAGCCGAGATAAGCCACACAAGACGCATGTTGGGTGGCTCACCGCATTATCAT
GCCACTCCATCAGGCCAGCCTGAGGAGAGATGAGGGCCGTGGACCAGGTGAATGCAAGCAAGGTGG
TCCGGAAGCACTCAGGGACAGACCGAACTGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR204497 protein sequence Red=Cloning site Green=Tags(s) MAATNLENQLHSAQKNLLFLQREHASTLKGLHAEIRRLQQHCTDLTYELTLKSFELTGDSSSRRTTELKRR CEELEAQLKAKEEENRKLLQEQKNAAIAVLENTVREREKKYLEELKVKSHKLSMLSGELEQRASTVAY LTSQLHAAKKKLLSSSGTSDASPAGSPALASYKPTPPKDKLPETPRRRMKKSLSAPLHPEFEEVYRFGAE SRKLLLRPVDAMPDPTPFLARESAEVQLKERPLVIPPIASDRSATGQHSPARDKPHKTHVGVAHRIHH ATPSQAQPEGEMRAVDQVNASKVVRKHSGTDRTV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_144819
ORF Size:	942 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_144819.2, NP_659068.1</u>
RefSeq Size:	2256 bp
RefSeq ORF:	945 bp
Locus ID:	215707
UniProt ID:	<u>Q8VDN4</u>
Cytogenetics:	5 63.85 cM
MW:	35.2 kDa