

Product datasheet for **MR204528L3V**

Aagab (NM_025857) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Aagab (NM_025857) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Aagab
Synonyms:	2310007F21Rik; p34
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025857
ORF Size:	948 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204528).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025857.2
RefSeq Size:	2600 bp
RefSeq ORF:	951 bp
Locus ID:	66939
UniProt ID:	Q8R2R3
Cytogenetics:	9 C
Gene Summary:	May be involved in endocytic recycling of growth factor receptors such as EGFR. [UniProtKB/Swiss-Prot Function]



[View online »](#)