

Product datasheet for MR204546

Ppp1r3c (NM_016854) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Ppp1r3c (NM_016854) Mouse Tagged ORF Clone
 Tag: Myc-DDK
 Symbol: Ppp1r3c
 Synonyms: Ppp1r5; PTG
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 ORF Nucleotide Sequence: >MR204546 representing NM_016854
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGGATCGCC

ATGAGCTGCACCAGAATGATCCATGTGCTAGATCCACGTCCTTTGACAAGTTCCGTCATGCCCGTGGACA
 TGGCCATGAGGATTTGCTTGGCTCATTACCACCTCTGAAGAGTTTCTGGGTCTTACAATGGTTTTCA
 ACGAAGAAATTTGTGAATAAATTGAAACCTTTGAAACCATGTCTCAGTGTCAAGCAGGAAGCCAAATCG
 CAGAGTGAGTGAAGAGCCACACAACCAAGCCAAGAAGCGGGTCTGTGTTGCGGACTCCAAGGGGCTGT
 CACTACTGCTATCCATGTCTTCTCCGACCTTCCAGAAGAACCAGCGTGGGACCTGCAGTTTGATCTCTT
 GGACCTTAACGATATCTCTCCAGCTTAAACTTCACGAGGAGAAAAATTTGGTTTTGATTTCCCCAG
 CCCTCAACCGACTACTTAAGTTTCCGGGACCGCTTTCAGAAGAAGCTTTGTCTGCCTCGAGAAGTGTCTT
 TGCAAGATCGGACGGTGACCGGGACAGTGAAGTGAAGAATGTGAGCTTTGAGAAGAAGTTTCAGGTCGG
 GATCACCTTTGACACCTGGAAAACCTACACAGATGTGGACTGTGTCTACATGAAGAATGTTTACAGCAGC
 TCAGACAGCGACACCTTCTCTTTGCCATCGACTTGCCCCGTGTCACTCCAAGTGAAGGAGAAAAATTGAGT
 TCTGCATTTCTTATCACGCTAATGGGAGGATCTTCTGGGACAACAATGAGGGTCAGAATTACAGAATTGT
 CCATGTGCAATGGAAAACCTGACGGAGTGCAGACTCAGGTGGCACCCAAAGACTGTGCATTCCAACAGGGG
 CCCCCTAAGACTGAGATAGAGCCACAGTCTTTGGCAGTCCAAGGCTTGCTAGCGGCCTCTTCCAGAGT
 GGCAGAGCTGGGGGAGAGTGGAGAACTTGACCTCCTATCGA

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR204546 representing NM_016854
 Red=Cloning site Green=Tags(s)

MSCTRMIVHLDPRLTSSVMPVDMAMRICLAHSPPLKSFLGPYNGFQRRNFVNKLKPLKPCLSVKQEAKS
 QSEWKSPhNQAKKRvVfADSKGLSLTAIHVfSDLPEEPAWDLQFDLLDLNDISSSLKLHEEKLVDFPQ
 PSTDYLsFRDRfQKNfVCLenCSLQDRtVTGTvKVKNVsfEKkVQVRITfDTWkTYTDVDCvYMKNVYSS
 SDSDTfSFaIDLPRvIPTeEKIEfCISyHANGRIfWdNNEGQNYRIVHVQWKPDGVQTVAPKDCAFQQG
 PPKTEIEPTVFGSPRLASGLfPEWQSWGRVENLTSYR

TRTRPLEQKLISEEDLAANDILDYKDDDDKv

Chromatograms: https://cdn.origene.com/chromatograms/mm9059_a08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_016854


ORF Size: 951 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_016854.2</u> , <u>NP_058550.1</u>
RefSeq Size:	2696 bp
RefSeq ORF:	954 bp
Locus ID:	53412
UniProt ID:	<u>Q7TMB3</u>
Cytogenetics:	19 C2
MW:	36.5 kDa
Gene Summary:	Acts as a glycogen-targeting subunit for PP1 and regulates its activity. Activates glycogen synthase, reduces glycogen phosphorylase activity and limits glycogen breakdown. Dramatically increases basal and insulin-stimulated glycogen synthesis upon overexpression in a variety of cell types.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR204546