

Product datasheet for **MR204562A1V**

Mouse Rplp0 (BC011106) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Rplp0 (BC011106) AAV Particle
Tag:	Myc-DDK
Symbol:	Rplp0
Synonyms:	36B4; Arbp; L10E
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR204562 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCCAGGGAAGACAGGGCGACCTGGAAGTCCAACACTTCTCAAGATCATCCAACCTTTGGATGATT
ATCCAAATGCTTCATTGTGGGAGCAGACAACGTGGGCTCCAAGCAGATGCAGCAGATCCGCATGTCCT
CCGAGGGAAGGCCGTGGTGCTGATGGGCAAGAACCATGATGCGCAAGGCCATCAGGGGCCACCTGGAG
AACACCCAGCTCTGGAGAACTGCTGCCTCACATCCGGGGGAACGTGGGCTTCGTGTTACCAAGGAGG
ACCTCACTGAGATTCGGGATATGCTGTTGGCCAATAAGGTGCCAGCTGCTGCTCGGGCTGGTGCCATCGC
CCCGTGTGAGGTCACTGTGCCAGCTCAGAACACTGGTCTGGGGCCGAGAAGACCTCCTTCTTCCAGGCT
TTGGGCATCACCACGAAAATCTCCAGAGGCACCATTTGAAATTCTGAGTGATGTGCAGCTGATAAAGACTG
GAGACAAGGTGGGAGCCAGCGAGGCCACACTGCTGAACATGCTGAACATCTCCCCCTTCTCCTTCGGGCT
GATCATCCAGCAGGTGTTTGACAACGGCAGCATTTATAACCTGAAGTGCTCGACATCACAGAGCAGGCC
CTGCACTCTCGCTTCTGGAGGGTGTCCGCAACGTGGCCAGTGTGTCTGCAGATCGGGTACCAACTG
TTGCCTCGGTGCCCACTCCATCATCAATGGGTACAAGCGCTCCTGGCATTGTCTGTGGAGACTGAGTA
CACCTTCCCACTTGCTGAAAAGGTCAAGGCCTTCTGGCTGATCCATCTGCATTTCGGGCTGCTGCCCT
GCAGCTGCTGCCACCACTGCTGCCCTGCGGTGCTGCAGCCCTGCCAAAGCTGAAGCAAAGGAAGAGT
CGGAGGAATCAGATGAGGATATGGGATTCGGTCTCTTCGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR204562 protein sequence Red=Cloning site Green=Tags(s) MPREDRATWKSNYFLKIIQLDDYPKCFIVGADNVGSKQMQQIRMSLRGKAVVLMGKNTMMRKAIRGHLE NNPALEKLLPHIRGNVGFVFTKEDLTEIRDMLLANKVPAAARAGAIAPCEVTVPAQNTGLGPEKTSFFQA LGITTKISRGTIIEILSDVQLIKTGDKVGASEATLLNMLNISPFSGLI IQQVFDNGSIYNPEVLDITEQA LHSRFLEGVRNVASVCLQIGYPTVASVPHSIINGYKRVLALSVETEYTFPLAEKVKAFLADPSAFAAAAP AAAATTAAPAAAAAPAKAEAKEESESEDEDMGFGFLFD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	BC011106
ORF Size:	951 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>BC011106</u> , <u>AAH11106</u>
RefSeq Size:	1138 bp
RefSeq ORF:	953 bp
Locus ID:	11837
Cytogenetics:	5 F
MW:	34.2 kDa