

## Product datasheet for **MR204562L3V**

### Arbp (BC011106) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Arbp (BC011106) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arbp
Synonyms:	36B4; Arbp; L10E
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC011106
ORF Size:	951 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204562).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC011106</a> , <a href="#">AAH11106</a>
RefSeq Size:	1138 bp
RefSeq ORF:	953 bp
Locus ID:	11837
Cytogenetics:	5 F
Gene Summary:	Ribosomal protein P0 is the functional equivalent of E.coli protein L10.[UniProtKB/Swiss-Prot Function]



[View online »](#)