

Product datasheet for **MR204566L3V**

Tfec (NM_031198) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tfec (NM_031198) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tfec
Synonyms:	BB107417; bHLHe34; Tcfec
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031198
ORF Size:	954 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204566).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031198.1 , NP_112475.1
RefSeq Size:	1779 bp
RefSeq ORF:	954 bp
Locus ID:	21426
UniProt ID:	Q9WTW4
Cytogenetics:	6 A2



[View online »](#)

Gene Summary:

Transcriptional regulator that acts as a repressor or an activator. Acts as a transcriptional transactivator on the proximal promoter region of the tartrate-resistant acid phosphatase (TRAP) E-box containing promoter. Collaborates with MITF in target gene activation. Acts as a transcriptional repressor on minimal promoter containing element F (that includes an E-box sequence) (By similarity). Binds to element F in an E-box sequence-specific manner (By similarity). Acts as a transcriptional repressor on minimal promoter containing mu E3 enhancer sequence (By similarity). Binds to mu E3 DNA sequence of the immunoglobulin heavy-chain gene enhancer (By similarity). Binds DNA in a homo- or heterodimeric form. [UniProtKB/Swiss-Prot Function]