

## Product datasheet for **MR204568L3V**

### Eci3 (NM\_026947) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Eci3 (NM_026947) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Eci3
Synonyms:	1810022C23Rik; Kseg2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026947
ORF Size:	954 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204568).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026947.3</a>
RefSeq Size:	1346 bp
RefSeq ORF:	954 bp
Locus ID:	69123
UniProt ID:	<a href="#">Q78JN3</a>
Cytogenetics:	13 A3.3
Gene Summary:	Catalyzes the isomerization of trans-3-nonenoyl-CoA into trans-2-nonenoyl-CoA (PubMed:24344334). May also have activity towards other enoyl-CoA species (Probable). [UniProtKB/Swiss-Prot Function]



[View online »](#)