

## Product datasheet for MR204569

### Nifk (NM\_026472) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Nifk (NM_026472) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nifk
Synonyms:	A1852665; C130020J04Rik; Mki67ip
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR204569 representing NM_026472 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

ATGGCTGGGTTAGCAGGCCCGCGAAGCCGAGTCTGGCATTGAACCCGCAGGAAGACTCGCAGTTTGAGA  
AGGCTCTGACGCAGATCCAAGGGCGCACGAAGAAGCCTCAGCAAAAGAAGAAAGAAAACTTAACCGTGG  
AGTAGTTTATCTGGGCCACCTACCTTCAACACTTTCTGAAAGCCACATCTATAACTACTGTGCCAGTTT  
GGCGACATTAGCAGATTCAGATTGTCTAGAAGTAAGCGGACTGGAACAGCAAAGGCTATGCCTTTGTGG  
AGTTTGAGTCTGAGGATGTTGCCAAGATAGTTGCAGAAACAATGGATAACTACCTTTTTGGTGAAGACT  
TCTATCATGTAAATTTATGCCACGAAAAAAGTCCATAAAGACCTTTTCAGCCAGAGGAATGCTCTGTTT  
CACCGGCCATCATTCCCAGCAGTGAAGCGCTATAATCGGAAACGAGGACATTTGCAGATGTTGAAGATGG  
AGTATAGGTTTAAAGAAGAAGGAAAAATTAAGCCAAAGAACTAGCTGCAAGGGGATTGATTATAGTTT  
TCCCTCACTGGTCTACCTAAGCCAAAGAACATTGCTGTGCTCATAGAGATTCTGAGGGTAACCCAGGTT  
TTACCTGACCAAAAAGAAGGGCTTTTCGGGTGAACCTAGGAGAAAGGAGAAGATGATGAAAGAGGACATTT  
CCAATAACATTCTAAGAAGCGTAAGAGGTCTAGGAGGAAGAAGTCGTCTGTGGACAGCCAGGGTCCCAC  
ACCAGTTTGTACACCGACGTTTTTGGAGAGGCGAAAATCGCAGGTGATGGAAGTCGGTGGCGATAAAGAT  
GATGAAATCATTCTCAAGCTGCCTGTGCCCCAGTGAAAGAAGACACAAAAAGACCCCAACATCTGCAA  
GCCCTGGGGAAAAAGACCGAGAAAAAGGAAGAGCAAGCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: [https://cdn.origene.com/chromatograms/mm9077\\_h07.zip](https://cdn.origene.com/chromatograms/mm9077_h07.zip)

Restriction Sites: SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_026472

**ORF Size:** 951 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_026472.4](#), [NP\\_080748.3](#)

**RefSeq Size:** 1894 bp


**RefSeq ORF:** 954 bp

**Locus ID:** 67949

**UniProt ID:** [Q91VE6](#)

Cytogenetics: 1 E2.3  
 MW: 36.7 kDa

**Product images:**



Circular map for MR204569