

Product datasheet for **MR204576A1V**

Mouse Ube2j1 (NM_019586) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Ube2j1 (NM_019586) AAV Particle
Tag:	Myc-DDK
Symbol:	Ube2j1
Synonyms:	0710008M05Rik; 1110030I22Rik; Ncube; NCUBE-1; Ncube1; Ubc6p
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR204576 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGACCCGCTACAACCTGAAGAGCCCGCCGTTAAGCGATTGATGAAAGAAGCTGCAGAATTGAAAG
ATCCAACAGATCATTACCACGCTCAGCCTCTGGAGGATAACCTTTTGAATGGCACTTCACAGTTAGAGG
GCCCCAGATTCGATTTTGATGGAGGAGTTTACCATGGACGCATAGTACTGCCACCAGAGTACCCCATG
AAACCACCAAGCATTATCCTGCTAACAGCTAATGGACGATTTGAAGTAGGCAAGAAGATCTGTTTAAGCA
TCTCAGGACATCATCTGAAACCTGGCAGCCTTCATGGAGCATAAGGACAGCCTTGCTGGCCATCATCGG
GTTTATGCCGACTAAAGGAGAGGGGCCATAGTTCTCTAGATTACACACCAGAGGAGAGAAGAGCGCTT
GCCAAAAAATCGCAAGATTTCTGTTGTGAAGGCTGTGGCTCTGCCATGAAGGATGTCTGTTGCCATTAA
AATCTGGAAGTGGTTCAAGCCAGGCTGACCAAGAAGCCAAAGAGCTGGCCAGACAGATTAGTTTTAAGGC
CGAAGTCAATTCATCTGGAAGACTATCGCTGAATCAGACTTAAACAGTGTTTCTCACTAAATGACTCA
CAAGACGATTTACCTACAACATTCCAGGGTGCTACAGCAAGCACATCGTATGGCGCCAGAATCCATCAG
GAGCACCTTGCCACAACCTACCCAGCCAGCTCCTAAGAATACCTCCATGAGCCCTCGGCAGCGCAGAGC
CCAGCAACAGAGTCAGAGAAGGCCATCCACTTCCCCGATGTGCTCCAAGGCCAGCCACCAAGAGCCAC
CACACAGAGCACGGCGGGTCCGCCATGCTGATCATCATCTGACTTTAGCATTGGCGGCTTTATATTCC
GACGAATATATCTGGCCAATGAGTACATATTTGACTTTGAGTTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR204576 protein sequence Red=Cloning site Green=Tags(s) METRYNLKSPAVKRLMKEAAELKDPTDHYHAQPLEDNLFWEHFTVRGPPDSDFDGGVYHGRIVLPPEYPM KPPSIILLTANGRFEVGKKICLSISGHPETWQPSWSIRTALLAIIGFMPTKGEGAIGSLDYTPERRAL AKKSQDFCCEGCGSAMKDVLLPLKSGSGSSQADQEAKELARQISFKAENVSSGKTIAESDLNQCFSLNDS QDDLPTTFQGATASTSYGAQNPSGAPLPQPTQPAPKNTSMSPRQRRRAQQSQRRPSTSPDVLQGQPPRAH HTEHGGSAMLIILLTLALAALIFRRIYLANEYIFDFEL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_019586
ORF Size:	954 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_019586.3 , NP_062532.2
RefSeq Size:	3522 bp
RefSeq ORF:	957 bp
Locus ID:	56228
UniProt ID:	Q9JJZ4
Cytogenetics:	4 A5
MW:	35 kDa