

## Product datasheet for **MR204580L4V**

### **Acy3 (NM\_027857) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Acy3 (NM_027857) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Acy3
Synonyms:	0610006H10Rik; AA3; AAIII; Acy-3; AW107362; HCBP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_027857
ORF Size:	957 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204580).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_027857.3</a>
RefSeq Size:	1505 bp
RefSeq ORF:	957 bp
Locus ID:	71670
Cytogenetics:	19 A



[View online »](#)

**Gene Summary:**

This gene encodes a member of the aminoacylase family of enzymes. This enzyme specifically deacetylates N-acetyl aromatic amino acids and mercapturic acids. Action of this enzyme on metabolites of the environmental contaminant trichloroethylene leads to the generation of toxic products that may lead to kidney failure. This protein has been found to bind to the hepatitis C virus core protein. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2014]