

## Product datasheet for **MR204588L3V**

### Slain1 (BC025223) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slain1 (BC025223) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slain1
Synonyms:	9630044O09Rik; AA675320; AW742596
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC025223
ORF Size:	954 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204588).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC025223</a> , <a href="#">AAH25223</a>
RefSeq Size:	1213 bp
RefSeq ORF:	956 bp
Locus ID:	105439
Cytogenetics:	14 E2.3
Gene Summary:	Microtubule plus-end tracking protein that might be involved in the regulation of cytoplasmic microtubule dynamics, microtubule organization and microtubule elongation. [UniProtKB/Swiss-Prot Function]


[View online »](#)