

Product datasheet for **MR204592**

Prps2 (BC024942) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Prps2 (BC024942) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Prps2
Synonyms:	MGC28388
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR204592 representing BC024942 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCTAACATCGTGCTCTTCAGCGGCAGTTCTCACCAGGACCTGTCGCAGCGCGTGGCCGATCGCCTGG
GCCTGGAGCTGGCAAGGTGGTCACCAAGAAGTTAGCAACCAGGAGACCAGTGTGAGATCGGTGAAAAG
TGTGAGAGGGGAAGACGTCTACATCATCCAGAGTGGCTGCGGGGAGATTAATGACAACCTGATGGAGCTG
CTCATCATGATCAATGCGTGCAAGATTGCGTCCTCATCCAGGGTGACAGCTGTGATACCTGCTTCCCCT
ACGCCCGCAAGATAAAAAGGACAAGAGTCGTGCTCCAATTTCTGCTAACTTGTGGCCAATATGCTGTC
AGTTGCTGGGGCTGATCACATCATTACCATGGACCTGCATGCCTCTCAGATCCAGGGTTTCTTTGATATT
CCTGTGGATAATTTATATGCAGAACCTGCAGTTCTTCAGTGGATTCCGGGAGAATCACTGAATGGAGAA
ACTGTATTATTGTTTCTCCAGATGCTGGAGGAGCCAAAAGGGTTACATCAATTGCAGACAGACTGAATGT
GGAATTTGCTTTGATTCACAAAGAGAGAAAAGAAAGCAAAATGAAGTGGACCGGATGGTTCTCGTGGGAGAT
GTGAAGGACCGTGTGGCCATCCTGGTAGACGACATGGCCGATACATGTGGTACCATCTGCCATGCTGCAG
ATAAGTTACTCTCAGCTGGTGCCACCAAGTTTATGCTATTCTTACCATGGAATCTTCTCTGGGCCAGC
TATTTCCAGAATAAACAGTGTGCCTTTGAAGCTGTTGTTGTACCAACACAATTCACAAGAGGACAAA
ATGAAACTGCTCCAAAATTCAGGTGATTGACATTTGATGATCCTGGCTGAAGCCATTGGAAGAACC
ACAATGGGAATCCGTGTCTATCTGTTCCAGCCATGTCCCACTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR204592 representing BC024942
 Red=Cloning site Green=Tags(s)

MPNIVLFGSSSHQDL SQRVADRLGLELGKVVTKKFSNQETSVEIGESVRGEDVYIIQSGCGEINDNLMEL
 LIMINACKIASSSRVTAVIPCFPYARQDKKDSRAPISAKLVANMLSVAGADHIIITMDLHASIQGFFDI
 PVDNL YAEPV LQWIRE NIT EWRNCII VSPDAGGAKRVT SIADRLNVEF ALIHKERKKANEVDRMVLVGD
 VKDRVAILVDDMADTCGTICHAADKLLSAGATKVYAILTHGIFSGPAISRINSAAFEAVVVTNTIPQEDK
 MKHCSKIQVIDISMILAEAIRRTHNGESVSYLFSHVPL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9062_h06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC024942

ORF Size: 954 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC024942.1](#)

RefSeq Size: 2860 bp

RefSeq ORF: 956 bp


Locus ID: 110639

Cytogenetics: X 78.32 cM

MW: 34.8 kDa

Gene Summary: Catalyzes the synthesis of phosphoribosylpyrophosphate (PRPP) that is essential for nucleotide synthesis.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR204592