

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_021554.2 , NP_067529.2
RefSeq Size:	1850 bp
RefSeq ORF:	957 bp
Locus ID:	59052
UniProt ID:	Q9EPL4
Cytogenetics:	7 F2
Gene Summary:	Protein-histidine N-methyltransferase that specifically catalyzes 1-methylhistidine (pro-methylhistidine) methylation of target proteins (PubMed:33563959). Mediates methylation of proteins with a His-x-His (HxH) motif (where 'x' is preferably a small amino acid) (PubMed:33563959). Catalyzes methylation of target proteins such as S100A9, NDUFB3, SLC39A5, SLC39A7, ARMC6 and DNAJB12; 1-methylhistidine modification may affect the binding of zinc and other metals to its target proteins (PubMed:33563959). Constitutes the main methyltransferase for the 1-methylhistidine modification in cell (PubMed:33563959). [UniProtKB/Swiss-Prot Function]

