

Product datasheet for **MR204595L3V**

Hnrpa3 (BC062198) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hnrpa3 (BC062198) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hnrpa3
Synonyms:	MGC37309, hnRNPA3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC062198
ORF Size:	954 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204595).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC062198 , AAH62198
RefSeq Size:	1781 bp
RefSeq ORF:	956 bp
Locus ID:	229279
Cytogenetics:	2 C3
Gene Summary:	Plays a role in cytoplasmic trafficking of RNA. Binds to the cis-acting response element, A2RE. May be involved in pre-mRNA splicing (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)