

Product datasheet for **MR204616L4V**

Wdr78 (BC034130) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Wdr78 (BC034130) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Wdr78
Synonyms:	MGC31722, MGC38950
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC034130
ORF Size:	960 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204616).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC034130 , AAH34130
RefSeq Size:	2560 bp
RefSeq ORF:	962 bp
Locus ID:	242584
Cytogenetics:	4 C6
Gene Summary:	Plays a critical role in the assembly of axonemal dynein complex, thereby playing a role in ciliary motility.[UniProtKB/Swiss-Prot Function]



[View online »](#)