

Product datasheet for **MR204705A1V**

Mouse Ccnh (NM_023243) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Ccnh (NM_023243) AAV Particle
Tag:	Myc-DDK
Symbol:	Ccnh
Synonyms:	6330408H09Rik; AI661354; AV102684; AW538719
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR204705 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTATCACAGCAGCAGCCAGAAACGGCACTGGACCTTCGCTAGCGAGGAGCAGCTGGCGCGCTCGGG
CGGACGCCAACCGCAAATTCAGTGCAAAGCTGTGGCTAACGGGAAGGTTCTTCCAAATGATCCAGTGTT
TCTGGAGCCTCATGAAGAATTGACACTTTGCAAGTACTATGAAAAAGATTATTGGAATTCTGTTCAAGT
TTTAAACCAGCTATGCCACGGTCTGTTGTGGGTACAGCTTGATGTATTTCAAGCGTTTTATCTTAATA
ACTCAGTAATGGAATATCACCTCGGATAATAATGCTTACTTGTGCATTTTGGCCTGCAAAGTAGATGA
ATTCAATGTGTCTAGTCCCCAGTTTGTGGAAATCTTCGAGAGAGTCCTCTGGACAGGAGAGGGCACTG
GAACAGATTTTGAATATGAACTACTACTTATACAACAACTTAATTTTACCTTATTGTCCACAATCCAT
ACAGACCATTGGAAGGCTTCCTCATCGATATAAAGACTCGGTACCCCATGTTGGAGAATCCGGAGATTCT
GAGGAAAACAGCTGATGATTTTCTAGTAGAATTGCATTGACAGATGCTTATCTTTATACACACCTTCA
CAGATTGCCCTGACCGCCATTTATCAAGTGCCTCTAGAGCTGGAATTACTATGGAAAGCTACTTATCAG
AGAGTCTAATGCTAAAAGAAAACAGAACTTGCCTGTACAGTTACTGGATATAATGAAAAGCATGAGAAA
TCTAGTAAAAAGTACGAGCCGCCAGGTCTGACGAGGTTGCTGTCTGAAACAGAAAGCTGGAGCGGTGT
CATTCTTCTGACCTTGCACTTAATGCGGTTACGAAGAAGAGGAAAGGCTATGAAGATGATGACTATGTGT
CAAAGAAACCCAAACAGGAAGAGGAGGAATGGACTGATGACGACCTGGTAGATTCTCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR204705 protein sequence Red=Cloning site Green=Tags(s) MYHSSSQKRHWTFASEEQLARLRADANRKFCKAVANGKVLNDPVFLEPHEELTLCKYYEKRLLEFCSV FKPAMPRSVVGTACMYFKRFYLNNSVMEYHPRIIMLTCAFLACKVDEFNVSSPQFVGNLRESPLGQERAL EQILEYELLLIQQLNFHLIVHNPYRPFEGFLIDIKTRYPMLENPEILRKTADDFLSRIALTDAYLLYTPS QIALTAILSSASRAGITMESYLSLMLKENRTCLSQLLDIMKSMRNLVKKYEPSPRDEVAVLKQKLERC HSSDLALNAVTKKRKGIEDDYVSKPKQEEEEWTDDDLVDSL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_023243
ORF Size:	969 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_023243.5</u>
RefSeq Size:	1976 bp
RefSeq ORF:	972 bp
Locus ID:	66671
UniProt ID:	<u>Q61458</u>
Cytogenetics:	13 C3
MW:	37.5 kDa