

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for MR204706L4V

Cd47 (NM_010581) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cd47 (NM_010581) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cd47
Synonyms:	9130415E20Rik; AA407862; Al848868; AW108519; B430305P08Rik; IAP; Itgp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010581
ORF Size:	972 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204706).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<u>NM 010581.3</u>
RefSeq Size:	1928 bp
RefSeq ORF:	975 bp
Locus ID:	16423
UniProt ID:	<u>Q61735</u>
Cytogenetics:	16 B5



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Cd47 (NM_010581) Mouse Tagged ORF Clone Lentiviral Particle – MR204706L4V

Gene Summary:Has a role in both cell adhesion by acting as an adhesion receptor for THBS1 on platelets, and
in the modulation of integrins. Plays an important role in memory formation and synaptic
plasticity in the hippocampus. Receptor for SIRPA, binding to which prevents maturation of
immature dendritic cells and inhibits cytokine production by mature dendritic cells.
Interaction with SIRPG mediates cell-cell adhesion, enhances superantigen-dependent T-cell-
mediated proliferation and costimulates T-cell activation. May play a role in membrane
transport and/or integrin dependent signal transduction. May prevent premature elimination
of red blood cells. May be involved in membrane permeability changes induced following
virus infection (By similarity).[UniProtKB/Swiss-Prot Function]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US