

## Product datasheet for **MR204725L3V**

### Spdef (NM\_013891) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Spdef (NM_013891) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Spdef
Synonyms:	PDEF; Pse
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013891
ORF Size:	978 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204725).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013891.2</a>
RefSeq Size:	1824 bp
RefSeq ORF:	978 bp
Locus ID:	30051
UniProt ID:	<a href="#">Q9WTP3</a>
Cytogenetics:	17 A3.3



[View online »](#)

**Gene Summary:**

May function as an androgen-independent transactivator of the prostate-specific antigen (PSA) promoter. Binds to 5'-GGAT-3' DNA sequences. May play a role in the regulation of the prostate gland and/or prostate cancer development. Acts as a transcriptional activator for SERPINB5 promoter (By similarity).[UniProtKB/Swiss-Prot Function]