

Product datasheet for **MR204726L2V**

Ccdc49 (BC079553) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ccdc49 (BC079553) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ccdc49
Synonyms:	1300013D05Rik; AA473636; Ccdc49; R75228
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	BC079553
ORF Size:	975 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204726).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC079553 , AAH79553
RefSeq Size:	5862 bp
RefSeq ORF:	977 bp
Locus ID:	67480
Cytogenetics:	11 D
Gene Summary:	Involved in pre-mRNA splicing as component of the spliceosome.[UniProtKB/Swiss-Prot Function]



[View online »](#)