

Product datasheet for **MR204739L1V**

Alx4 (BC049786) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Alx4 (BC049786) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Alx4
Synonyms:	Ist
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	BC049786
ORF Size:	975 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204739).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC049786
RefSeq Size:	7133 bp
RefSeq ORF:	977 bp
Locus ID:	11695
Cytogenetics:	2 51.62 cM
Gene Summary:	Transcription factor involved in skull and limb development.[UniProtKB/Swiss-Prot Function]



[View online »](#)