

Product datasheet for **MR204744L4V**

Adora1 (BC079624) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adora1 (BC079624) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adora1
Synonyms:	A1AR, A1R
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC079624
ORF Size:	978 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204744).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC079624 , AAH79624
RefSeq Size:	4607 bp
RefSeq ORF:	980 bp
Locus ID:	11539
Cytogenetics:	1 E4
Gene Summary:	Receptor for adenosine. The activity of this receptor is mediated by G proteins which inhibit adenylyl cyclase.[UniProtKB/Swiss-Prot Function]



[View online »](#)