

Product datasheet for **MR204777L3V**

Apof (NM_133997) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Apof
Synonyms:	AI255964
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PSI00092)
Tag:	Myc-DDK
ACCN:	NM_133997
ORF Size:	978 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR204777).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_133997.2 , NP_598758.1
RefSeq Size:	1930 bp
RefSeq ORF:	981 bp
Locus ID:	103161
UniProt ID:	Q91V80
Cytogenetics:	10 76.5 cM



Gene Summary:

Minor apolipoprotein that associates with LDL. Inhibits cholesteryl ester transfer protein (CETP) activity and appears to be an important regulator of cholesterol transport. Also associates to a lesser degree with VLDL, Apo-AI and Apo-AII.[UniProtKB/Swiss-Prot Function]