

## Product datasheet for **MR204778L4V**

### Cd200r1 (NM\_021325) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cd200r1 (NM_021325) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cd200r1
Synonyms:	CD200R; Mox2r; OX2R
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021325
ORF Size:	981 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204778).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021325.2</a>
RefSeq Size:	1891 bp
RefSeq ORF:	981 bp
Locus ID:	57781
UniProt ID:	<a href="#">Q9ES57</a>
Cytogenetics:	16 B4



[View online »](#)

**Gene Summary:**

Inhibitory receptor for the CD200/OX2 cell surface glycoprotein. Limits inflammation by inhibiting the expression of proinflammatory molecules including TNF-alpha, interferons, and inducible nitric oxide synthase (iNOS) in response to selected stimuli.[UniProtKB/Swiss-Prot Function]