

Product datasheet for **MR204794L3V**

B4galt7 (NM_146045) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	B4galt7 (NM_146045) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	B4galt7
Synonyms:	MGC28643
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_146045
ORF Size:	984 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204794).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_146045.1 , NP_666157.1
RefSeq Size:	1636 bp
RefSeq ORF:	984 bp
Locus ID:	218271
UniProt ID:	Q8R087
Cytogenetics:	13 B1
Gene Summary:	Required for the biosynthesis of the tetrasaccharide linkage region of proteoglycans, especially for small proteoglycans in skin fibroblasts.[UniProtKB/Swiss-Prot Function]



[View online »](#)