

Product datasheet for **MR204798L4V**

Acp5 (NM_007388) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Acp5 (NM_007388) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Acp5
Synonyms:	TRACP; TRAP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007388
ORF Size:	981 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR204798).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007388.3 , NP_031414.1
RefSeq Size:	1367 bp
RefSeq ORF:	984 bp
Locus ID:	11433
UniProt ID:	Q05117
Cytogenetics:	9 8.38 cM
Gene Summary:	May play a role in the process of bone resorption. The osteoclastic trap acts on nucleotide tri- and diphosphates with higher affinity, compared with other substrates.[UniProtKB/Swiss-Prot Function]


[View online »](#)