

Product datasheet for **MR204806A1V**

Mouse Trib3 (BC012955) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Trib3 (BC012955) AAV Particle
Tag:	Myc-DDK
Symbol:	Trib3
Synonyms:	Trb3, Nipk, SINK, SKIP3, TRB-3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR204806 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCGAGCTACACCTCTGGCTGCTTCTGCTGATGTTTCCTGCAGGAAGAAACCGTTGGAGTTTGATGACA
ATATTGATGCCAAGTGTCCAGTCCTAAAACGAGTGAGAGATGAGCCTGAGCCCGGACCACTCCCCAGCCT
GCTGCCCCCAGCCACCTCCCGCCTCAGACTTGTCACCTGCTGTGGCCCTGCAACTCGACTGGGGCCT
TATATCCTTTTGAACGAGAGCAAGGCAGCTGCAGCTATCGAGCCCTGCACTGCCCCACAGGCACAGAGT
ACACCTGCAAGGTGTACCCTGCCAGCGAGGCCAGGCGGTGCTGGCACCTTATGCCCGGCTGCCTACCCA
CCAGCATGTGGCCCGTCCACAGAGGTCTGCTGGGCTCTCGGCTCCTTTACATCTTTTTCACGAAGACC
CATGGGACTTGACAGCCTGGTGCGCAGCCGCGCGGTATCCAGAGCCCGAGGCTGCCGGGCTCTTCC
GGCAGATGGCTAGTGCCGTGGCACACTGCCACAAGCACGGCTTGTCTTGCGCGACCTCAAGCTGCGTCG
CTTTGTCTTCAGCAACTGTGAGAGGACGAAGCTGGTGTGGAGAACCTGGAAGATGCCTGCGTGATGACT
GGATCAGATGACTCTCTGTGGGACAAGCATGCGTGCCCTGCCTACGTGGGACCAGAGATACTCAGTCCC
GGCCATCCTACTCTGAGCCAGTCTGCTCTTTGGCAAGATCCGTAGAGGGACCTTTGCCCTGCCTGAGGG
CCTATCAGCCCCAGCCGCTGTCTGATCCGCTGTCTCCTCCGAAGAACCTTCAGAGCGACTTGTGGCC
CTGGGAATCCTCTTGCATCCCTGGTTGAGAGAGGATCACGGCCGAGTCTCTCCTCCACAGTCTGACCGAA
GGGAGATGGACCAGGTGGTCCCAGATGGGCCACAGCTGGAGGAGGCTGAGGAAGGGGAGGTGGGGCTGTA
CGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR204806 protein sequence Red=Cloning site Green=Tags(s) MRATPLAASADVSCRKKPLEFDDNIDAKCPVLKRVRDEPEPGPLPSLLPPSPPPASDLSPAVAPATRLGP YILLEREQGSCSYRALHCPTGTEYTCVVYPASEAQAVLAPYARLPTHQHVARPTEVLLGSRLLYIFFTKT HGDLHSLVRSRRGIPEPEAAGLFRQMASAVAHCHKHGLVLRDLKLRRFVFSNCERTKLVLNLEDACVMT GSDDSLWDKHACPAYVGPEILSSRPSYSEPVLLFGKIRRGTFALPEGLSAPARCLIRCLLRKEPSERLVA LGILLHPWLREDHGRVSPPQSDRREMDQVVPDGPQLEEAEEGEVGLYG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	BC012955
ORF Size:	984 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>BC012955</u> , <u>AAH12955</u>
RefSeq Size:	1969 bp
RefSeq ORF:	986 bp
Locus ID:	228775
Cytogenetics:	2 74.83 cM
MW:	36.2 kDa