

Product datasheet for **MR204814**

Dmwd (BC031193) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dmwd (BC031193) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Dmwd
Synonyms:	DMR-N9, MGC37679, 59
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR204814 representing BC031193 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTTCTGCGTGGGCTCATGAAGAGCTACTTTGGGGTCTGCTGTGTGTGCTGGAGCCCTGATGGCC
GCTACGTGGTGACAGGAGGGGAAGATGACTTAGTCACTGTGTGGTCCTTCACGGAGGGTCGTGTGGTGGC
TCGAGGCCATGGCCACAAGTCTGGGTCAACGCTGTAGCCTTTGATCCCTACACCACTCGCGCAGAGGAG
GCGGCGTCAGCCAGTGTGATGGAGATCCAGTGGTGAAGAGGAGGAGCCTGAAGTACCAGCTCAGACA
CAGGAGCCCCAGTGTCCCACTGCCAAGGCTGGCTCCATCAGTACCGCTTTGGCTCAGCTGGCCAGGA
CACACAGTTCTGCCTGTGGGACCTCACAGAAGATGTGCTCTCCCTCATCCGTCTCTGGCCCGTACCCGC
ACCTTCCGGGCACACCTGGTGCCACCCACCAGCTTCTGGTAGTTCTCGGGCCGAGAGACAGGTGCAG
GCCCCCTGCCCGCTCCCTGTCTCGTTCCAACAGTCTCCACACCCAGCTGGTGGTGGCAAGGCTGGTGG
GCCTAGTGCATCGATGGAGCCTGGCATACCATTCAGCATTGGCCGCTTTGCCCACTGACCTGCAGGAG
CGGCGGGACCGGGGAGCTGAGAAGGAACACAAACGCTACCATAGCCTGGGAAACATCAGCCGCGTGGCA
GTGGGGCAATAGCAGCAATGACAAGCTCAGTGGTCTGCCCCCGAAGCCGATTGGACCCAGCTAAGGT
GCTGGGCACGGCACTGTGCCCTCGGATCCATGAGGTGCCACTGCTGGAGCCTCTCGTGTGCAAGAAGATT
GCTCAGGAACGCTGACCGTGCTACTGTTCTGGAGGATTGTATCATCACTGCCTGCCAAGAGGGCCTCA
TCTGCACCTGGGCCCGGCCAGGCAAGGCGGCATCTCTCCCAACCAGGCAGCTCCCCAGTGGCACTGT
GGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR204814 representing BC031193
 Red=Cloning site Green=Tags(s)

MLLRGLMKSYPFGLLCVCWSPDGRYVVTGGEDDLVTWVSFTEGRVVARGHGHKSWVNAVAFDPYTTAAEE
 AASASADGDPGSGEEPEVTSSDTGAPVSPLPKAGSITYRFGSAGQDTQFCLWDLTEDVLSHPHSLARTR
 TLPGTGATPPASGSSRAGETGAGPLPRSLSRNSLPHPAGGGKAGGPSASMEPGIPFSIGRFATLTLQE
 RRDRGAKEHKRYHSLGNISRGSGGNSNDKLGPAPRSRLDPAKVLGTALCPRIHEVPLLEPLVCKKI
 AQERLTVLLFLEDCIITACQEGLICTWARPAGKAGISSQPGSSPSGTVV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC031193

ORF Size: 984 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC031193.1](#)

RefSeq Size: 1629 bp


RefSeq ORF: 986 bp

Locus ID: 13401

Cytogenetics: 7 9.46 cM

MW: 34.3 kDa

Product images:



Circular map for MR204814