

Product datasheet for **MR204815A1V**

Mouse Clec11a (NM_009131) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Clec11a (NM_009131) AAV Particle
Tag: Myc-DDK
Symbol: Clec11a
Synonyms: AW457320; Clecsf3; Scgf
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR204815 representing NM_009131
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCAGGCAGCCTGGCTCTTGGGGGCCCTAGTGGTCCCTCAGCTTTTGAGTTTTGGTCATGGAGCCCGAG
GTCCTGGGAGGGAGTGGGAGGGAGGCTGGGGAGGTGCCCTGGAGGAGGAGAGAGAGCGGGAGTCACAGAT
GTTGAAGAATCTCCAGGAGGCCCTAGGGCTGCCCACTGGGGTGGGAAATGAGGATAATCTTGCTGAAAAC
CCTGAAGACAAAGAGGTCTGGGAGACCACAGAGACTCAAGGGGAAGAAGAGGAAGAGGAAATCACCACAG
CACCTTCTTCTAGTCCCAACCTTTCCCGAGCCCTTCTCCACACCCAGAGGACACTGTCACCTTACATCTT
GGGCCGCTTGGCCAGCCTCGATGCAGGCCTACACCAATTGCACGTCCGTCTGCACGTTTTGGACACCCGT
GTGTTGAGCTGACCCAGGGGCTGCGGCAGCTGCGGGATGCTGCGAGTGACACCCGCGACTCAGTGCAAG
CCCTGAAGGAGGTCCAGGACCGTCTGAGCAGGAGCACGGCCGCTTGGAGGGCTGCCTGAAGGGCCTGCG
CCTTGGCCACAAGTGCTTCTGCTCTCGCGAGACTTCGAGACCCAGGCGCGCGCAGGCGCGGTGCAAG
GCGCGAGGTGGGAGCTTAGCACAGCCTGCGGACCGCCAGCAAATGGATGCGCTAAGCCGGTACTTACGCG
CCGCTCTCGCCCCCTACAACCTGGCCGGTGTGGCTGGGAGTGCACGATCGGCGCTCCGAGGGGCTCTACCT
TTTCGAGAACGGCCAGCGGTGTCTTTCTCGCCTGGCACCAGCGCATTACGCTGGAGTCCGGCGCCAG
CCTAGTGCGGCAACACATCCACTCAGCCCGGATCAGCCAATGGCGGCGTCTTGAGAACTGCGTGGCCC
AGGCCTCAGACGACGGTTCTTGGTGGGACCATGACTGTGAGCGGCGCCTCTACTTCGTCTGCGAGTTCCC
CTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR204815 representing NM_009131 Red=Cloning site Green=Tags(s) MQAAWLLGALVVPQLLSFGHGARGPGREWEGGWGGALEEEERERESQMLKNLQEALGLPTGVGNEDNLAEN PEDKEVWETTETQGEIEEITAPSSSPNPFSPSPTPEDTVTYILGRLASLDAGLHQLHVRLHVLDT VVELTQGLRQLRDAASDTRDSVQALKEVQDRAEQEHGRLEGCLKGLRLGHKCFLLSRDFETQAAAQARCK ARGGSLAQPADRQQMDALSRYLRAALAPYNWPVWLGVDHRRSEGLYLFENGQRVSFFAWHRAFSLSEGAQ PSAATHPLSPDQPNGGVLENCVAQASDDGSWWDHDCERRLYFVCEFPF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_009131
ORF Size:	984 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_009131.3</u> , <u>NP_033157.1</u>
RefSeq Size:	1970 bp
RefSeq ORF:	987 bp
Locus ID:	20256
UniProt ID:	<u>O88200</u>
Cytogenetics:	7 B3
MW:	36.5 kDa