

Product datasheet for **MR204829**

BC024479 (BC024479) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BC024479 (BC024479) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	BC024479
Synonyms:	MGC37388
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR204829 representing BC024479 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGACTATTCACAGGAGGACTGGGAAACCACAGTCAGGAGCTCCATGGCTATCCAACAGATCAGG
AATTGGATGAAATGCCAGTTTCAAAGAGAACGTTAAAAATAAACAGAGTCTTCTGAAGAGGCACAGAA
GCGAGATACCATGCAGAAATTTGTACAGATTTTGGAGTCGGTGCAGTTGAAATGGGAAGTGTTCAGAGT
TGGACAGACTTCTCAAGGCTCCATCTTCCAATAAATTGGCCATTTTGAATTGGCTATAACACCCGCT
GGAAAGAGGATATTCGTTACCATTATGCTGAGATCAGCTCACAGGTGCCCTTGGCAAGCGTCTTCGCGA
GTATTTCAACTCTGAGAAAGCCTGAGGGACGGATCATTATGACCCGGGTGCAAAAAATGAAGTGGAAAAAT
GTTTACTACAAGTTTTTAGAGATCACCATCAGTGAAGCGAGGTGCCTAGAGCTGCACATGGAGATTGACT
GGATACCCATTGCCCCTCGAAGCCAACCGGGGCAATGTTGTTTCAGTATTTATTGCCAGGAGGCATTCC
TAAAGCCAGGCCTTTATGCCATTGGCTATGAAGAATGTATCGAGAGGCCTCTCTCACCAGATGTGGAG
CGCCATGCTCTGGACCCAGGGAAAGAGGGCCGAGTCGACCTGGAGACCCTCTCAGCACAAGCCTCATTAC
AGGTGGAAGTGAACCCACCCGGATCATCTACTGCTACCTCGGTATTGCTGAGGTGAGGACTCTGCAGCA
GTGCTTGTTTTACATTTCCAAGCAAATGCCAAACCTTCAGTAAAGAGTGGGTTGGTATTAATGGGTTT
TTGTCTCAGAACTGTATTGTGGATCCTGGAGTTTCCCCAAGTCCATCTACATCAAATTTGTAGAAGTAG
AGAGGGATTTTCTTCTGCCGGCTCATTAGTTGAGTGTCTGAAAAAGCCATTGGATACCCCTTAAAT
TAACAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR204829 representing BC024479
Red=Cloning site Green=Tags(s)



MEDYSQEDWGNHSQELHGYPDQELDEMPVSKRTLKIKQESSEEAQKRDTMQNIVQILESQVQLKWELFQS
WTDLSRLHLSNKLAIIFGIGYNTRWKEDIRYHYAEISSQVPLGKRLREYFNSEKPEGRIIMTRVQKMNWKN
VYKFLFLEITISEARCLELHMEIDWIPAHSKPTGGNVVQYLLPGGIPKSPGLYAIGYEECIERPLSPDVE
RHALDPGKEGRVDLETLSAQASLQVEVEPTRIIYCYLGIAEVRTLQQCLFLHFQANAKTFKEWVGINGF
LSQNCIVDPGVSPKSIYIKFVEVERDFLSAGSLVECLEKAIGYPLKFNN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: BC024479

ORF Size: 987 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC024479.1](#)

RefSeq Size: 1540 bp


RefSeq ORF: 989 bp

Locus ID: 235184

Cytogenetics: 9 A4

MW: 56.4 kDa

Product images:



Circular map for MR204829